




Infracrveni beskontaktni termometar za slepoočnicu HA500



Kutija za odlaganje



- rastojanje prilikom merenja je 5 cm
- brzo merenje u roku od jedne sekunde
- temperatura slepoočnice
- temperatura predmeta
- alarm povišene temperature
- memorija 9 merenja
- pozadinsko osvetljenje
- automatsko isključivanje
- indikator slabe baterije



Klinički tačan

ASTM E1965-98:2009

Ovaj termometar za slepoočnicu je klinički ispitan u nekoliko velikih nastavnih bolnica na osnovu protokola ASTM E1965-98: 2009 regulatorni standard, koji pokriva dovoljno povišene telesne temperature i normalne telesne temperature ispitanika sa zadovoljavajućom kliničkom ponovljivošću i tačnosti merenja u poređenju sa očitavanjem merenja oralne temperature.

SALVUS
INTERNATIONAL

Salvus International d.o.o.
Dobanovačka 56 | 11080 Zemun | Beograd
www.salvushealth.com | info@salvushealth.com