



# NI60 Qutie

Super mini kompresorski inhalator

**+** dizajnirano za decu



Tehnologija podesivog ventila

**Patentirani podesivi ventil s lakoćom daje lekove različitih stepena viskoznosti u skladu sa potrebama i stanjem svakog pacijenta.** Naša tehnologija podesivog ventila omogućuje korisnicima podešavanje različitih intenziteta raspršivanja i inhalacije u rasponu od 0,08 do 0,50 ml/min uz konstantnu veličinu čestica. Intenzitet inhalacije se podešava veoma jednostavno, bez zamene delova.




Zatvoreno



Potpuno otvoreno

- MMAD  $\leq 2,6 \mu\text{m}$ ;  
Udeo finih čestica: 75-80%
- Konzistentna veličina finih čestica za efikasan tretman respiratornog sistema
- Posuda za lek sa podesivim ventilom
- Nizak nivo buke
- Idealno za decu, pogodno i za odrasle
- Kompaktan uređaj male težine
- Nastavak za usta i maska za decu
- Adapter od 12 V



Efikasan tretman



Posuda za lek



Maska za decu



Nizak nivo buke



Nastavak za usta

Patentni brojevi :  
CN: ZL 2012 2 0443597.3  
ZL 2014 3 0053957.3  
Ostali patenti na čekanju

